

COSMEDENT®

Spreading Smiles™



CLASS III RESTORATION

Dr. Buddy Mopper



Class III



330 Carbide bur to begin prep



FlexiDiamond Strip used to clean interproximal contact



6862-012 diamond used to create long bevel



Notice the uniform 2mm long bevel



Mylar strip in place while applying Cosmedent's EnamelEtch, well past the long bevel using the #1 brush




Bonding adhesive in place (Cosmedent Complete)



Place Renamel Microhybrid A3 to increase value and opaque darkened area.



Sculpting of Renamel Microhybrid with Cosmedent's 8AL to the leading edge of the long bevel. Notice blending of color as well as increase in opacity



Renamel Microfill A3 rolled into a convenience form



Sculpting of Renamel Microfill layer using the Cosmedent 8AL instrument



Contouring of restoration using the Et9 carbide bur moving from composite to tooth



Use the Coarse FlexiDisc to blend the microfill margin into the tooth



Use the FlexiDisc 3/8th disc medium to finish at cervical



Use the FlexiDisc Medium 5/8th to contour, starting from the proximal and flexing to keep the proximal contour and create perfect embrasures



Use the FlexiDisc 5/8th Fine



Retraction of gingival tissue with the 8392-016 to protect tissue and reduce excess flash



Using the 8A plastic instrument to torque the teeth and gently separate the FlexiDimondStrip



Use the Fine/Super Fine FlexiStrip



Final polish with Enamelize and FlexiBuff disc.



Final Restoration